

Datasheet for S000-03**Streptavidin Peroxidase Conjugated****Overview**

Description:	Streptavidin Peroxidase Conjugated - S000-03
Item No.:	S000-03
Size:	1 mg
Applications:	Dot Blot, WB, ELISA, EM, IF, IHC, Other

Product Details

Background:	Streptavidin is isolated from bacteria, <i>Streptomyces avidinii</i> , and has an exceptionally high binding affinity for B7 (biotin). Rockland offers streptavidin in unconjugated and conjugated forms for common immunoassays including ELISA, western blotting, immunohistochemistry. Streptavidin is a tetrameric protein capable of binding 4 biotin groups to each molecule of streptavidin. While streptavidin has identical binding properties as avidin, it lacks the glycoprotein portion of the molecule and therefore shows less non-specific binding. Streptavidin is a slightly smaller molecule with a molecular weight of approximately 53.6 kDa. The sequence of avidin only shows 30% homology with streptavidin, and anti-avidin and anti-streptavidin antibodies are not immunologically cross reactive. Horseradish Peroxidase (HRP) is an enzyme that utilize organic peroxide compounds as electron donors. Naturally provides protection for plants against pathogens, but can be utilized in molecular biology to convert various substrates to detectable compounds (such as in Western Blotting and ELISAs).
Synonyms:	HRP-SA, Horseradish Peroxidase conjugated S avidin, Streptavidin HRP, Streptavidin conjugated to horseradish peroxidase (HRP), HRP-linked Streptavidin
Conjugate:	Peroxidase (HRP)

Target Details

Purity/Specificity:	Streptavidin-HRP was prepared from chromatographically purified streptavidin. Streptavidin Peroxidase conjugate was assayed by immunoelectrophoresis resulted in a single precipitin arc against anti-Peroxidase and anti-Streptavidin.
Relevant Links:	• S000-03 SDS • UniProtKB - P22629 • NCBI - CAA00084.1

Application Details

Tested Applications:	Dot Blot, WB
Suggested Applications:	ELISA, EM, IF, IHC, Other (Based on references)
Application Note:	Streptavidin Peroxidase has been tested by dot blot and western blot and is a useful detection reagent for primary antibodies conjugated to biotin. Streptavidin Peroxidase can be utilized in Immunohistochemistry, Immunofluorescence, immuno-EM, Western Blotting, and ELISA experiment formats in combination with the proper substrate (TMB-1000 or FEMTOMAX-110).
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.
ELISA:	1:20,000 - 1:200,000
IF:	User Optimized
IHC:	1:1,000 - 1:5,000
WB:	1:10,000 - 1:40,000


Formulation

Physical State:	Lyophilized
Concentration:	1.0 mg/mL by UV absorbance at 280 nm
Buffer:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2
Preservative:	0.01% (w/v) Gentamicin Sulfate. Do NOT add Sodium Azide!
Stabilizer:	10 mg/mL Bovine Serum Albumin (BSA) - Immunoglobulin and Protease free
Reconstitution Volume:	1.0 mL
Reconstitution Buffer:	Restore with deionized water (or equivalent)

Shipping & Handling


Shipping Condition:	Ambient
Storage Condition:	Store vial at 4° C prior to restoration. For extended storage aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. Streptavidin Peroxidase conjugated is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use.
Expiration:	Expiration date is one (1) year from date of receipt.

Images



Bottle

Streptavidin Peroxidase Conjugated



Dot Blot


Dot Blot Results of Streptavidin Peroxidase Conjugate.

Row A: BSA-Biotin Conjugated.

Row B: BSA.

Sample dilutions: 1- 100ng, 2- 33.33ng, 3- 11.11ng, 4-3.7ng, 5- 1.23ng.

Streptavidin Peroxidase Conjugated at 1.0µg/mL for 30mins at RT.



Western Blot

Western Blot Results using Streptavidin Peroxidase Conjugate and Goat Anti-GST Biotin Conjugate Antibodies.

Lane 1: Opal prestained molecular weight ladder (p/n MB-210-0500).


Lane 2: GST (p/n 000-001-200) [0.05 µg].

Primary Antibody: Goat Anti-GST Biotin Conjugate (p/n 600-106-200) at 1.0µg/mL overnight at 4°C.

Secondary Antibody; Streptavidin Peroxidase Conjugate at 1:40,000 for 30mins at RT.

Block: Blocking Buffer for Fluorescent Western Blotting (p/n MB-070) for 30mins at RT.


Exp: 5 sec.



Immunocytochemistry

Immuno-electron microscopy using biotinylated anti-rabbit and streptavidin-HRP.


TREM2^{-/-} mice show more severe axonal pathology after CPZ. (c) EM images of WT and TREM2^{-/-} at 12 weeks of CPZ treatment. Black arrows indicate dystrophic autophagocytic axons and asterisks indicate Iba1+ immunolabeled microglia. Fig. 2. PMID: 25631124.



Immunocytochemistry

Immuno-electron microscopy using biotinylated anti-rabbit and streptavidin-HRP.

Defect in myelin degradation in TREM2^{-/-} microglia. (c) Immuno-EM images of TREM2^{-/-} and WT microglia stained with Iba1 in the CC at 4 and 12 weeks on CPZ treatment. Images on the left in WT and TREM2^{-/-} panels at week 4 and 12 (3,000× magnification) depict Iba + microglial cells (asterisks). A higher magnification (15,000×) for the boxed area is shown on the right of each image. Black arrows indicate phagosomes containing myelin debris. White arrows indicate pi granules. Fig. 7. PMID: 25631124.




Western Blot

Western Blot of Goat anti-Glycerol Kinase Antibody Biotin Conjugated using Streptavidin HRP. Lane 1: Glycerol Kinase. Load: 50 ng per lane. Primary antibody: Glycerol Kinase Antibody Biotin Conjugated at 1:1000 overnight at 4°C. Secondary antibody: HRP Streptavidin (p/n S000-03) secondary antibody at 1:40,000 for 30 min at RT. Block: MB-070 for 30 min at RT. Predicted/Observed size: 55 kDa, 55 kDa for Glycerol Kinase.



Dot Blot

Dot Blot of Human IgG F(c) Fragment Biotin Conjugated using Streptavidin HRP. Human IgG F(c) Biotin Conjugate (1) 100ng, (2) 33.33ng, (3) 11.11ng, (4) 3.70ng, (5) 1.23ng. Primary Antibody: none. Secondary Antibody: Streptavidin HRP (p/n S000-03) at 1:40,000 for 30 mins at RT. Block: BlockOut buffer (p/n MB-073) at RT for 30 mins. Exposure: 1 sec.




ELISA

ELISA Results of Human IgG Whole Molecule Biotin Conjugated using Streptavidin-HRP. Each well was coated in duplicate with 1.0 µg of Human IgG Whole Molecule Biotin Conjugate. The working dilution is 82,800. The starting dilution of antibody was 5µg/ml and the X-axis represents the Log10 of a 3-fold dilution. This titration is a 4-parameter curve fit where the IC50 is defined as the titer of the antibody. Assay performed using Streptavidin-HRP (p/n S000-03) and TMB substrate (p/n TMBE-1000).



Immunohistochemistry

Percentage of cells showing cell-cycle activation as informed by Ki67 immunostaining (A) and cell death as identified by TUNEL assay (B) in the kidney in response to pregnancy in mice. Representative image of stained kidney from non-pregnant and pregnant female mice. Data are presented as mean ± SEM (n = 5/group). Asterisks represent significant differences between non-pregnant and pregnant mice as determined by Student's t-Test (* p < 0.05). Slides incubated with goat anti-rabbit secondary antibody (1:1000) and streptavidin-horseradish peroxidase (1:500, p/n S000-03). Images with the labels i and ii depict high magnification of the selected area. Arrow heads indicate positive DAB staining. Scale bar is 50 µm. Figure 3. PMID: 35682969.



Immunohistochemistry

The expression of p110 α protein by the mouse placenta on day 19 of pregnancy. Representative stained section shown with negative control shown in the figure inset. For localization of p110 α , placental sections were washed with PBS to remove OCT and underwent antigen retrieval with citrate buffer before immunolabelling against p110 α . Sections were treated with 0.5% Triton X-100 before immunolabelling. Bound antibody was detected using biotinylated goat anti-rabbit IgG followed by streptavidin-conjugated horseradish peroxidase (p/n S000-03) and 3,3'-diaminobenzidine (DAB). Sections were lightly counterstained with hematoxylin and mounted in DPX. Supplement Fig 1. PMID: 31241463.



Western Blot

Western Blot of Peroxidase Conjugated Streptavidin. Lane 1: Human IL-7. Load: 50 ng per lane. Primary antibody: Human IL-7 Biotin Conjugated antibody at 1:1,000 for overnight at 4°C. Secondary antibody: Peroxidase Conjugated Streptavidin at 1:40,000 for 30 min at RT. Block: 5% BLOTTO 30 min at RT. Predicted/Observed size: 17 kDa, 17 kDa for Human IL-7. Other band(s): none.

References

- Fial I et al. Characterizing Antibodies Targeting Antisense Oligonucleotide Phosphorothioate and 2'- O-Methoxyethyl Modifications for Intracellular Trafficking and Biodistribution Studies. *Nucleic Acid Ther.* (2025)
- Ruz-Maldonado I et al. Heterogeneity of hepatocyte dynamics restores liver architecture after chemical, physical or viral damage. *Nat Commun.* (2024)
- Sun J et al. Fatty acid binding protein 5 suppression attenuates obesity-induced hepatocellular carcinoma by promoting ferroptosis and intratumoral immune rewiring. *Nat Metab.* (2024)
- Lopez-Tello J et al. Fetal manipulation of maternal metabolism is a critical function of the imprinted Igf2 gene. *Cell Metab.* (2023)
- Mishra KP et al. Development of ELISA-Based Assay for Detection of SARS-CoV-2 Neutralizing Antibody. *Viral Immunol.* (2023)

- Ghezzi L et al. Schwann Cell Remyelination in the Multiple Sclerosis Central Nervous System. *Lab Invest.* (2023)
- Reuss S et al. Synaptoporin and parathyroid hormone 2 as markers of multimodal inputs to the auditory brainstem. *J Chem Neuroanat.* (2023)
- Lean SC et al. Obesogenic diet in mice compromises maternal metabolic physiology and lactation ability leading to reductions in neonatal viability. *Acta Physiol (Oxf).* (2022)
- Napso T et al. Diet-induced maternal obesity impacts feto-placental growth and induces sex-specific alterations in placental morphology, mitochondrial bioenergetics, dynamics, lipid metabolism and oxidative stress in mice. *Acta Physiol (Oxf).* (2022)
- Aykroyd BRL et al. Loss of imprinting of the Igf2-H19 ICR1 enhances placental endocrine capacity via sex-specific alterations in signalling pathways in the mouse. *Development.* (2022)
- Sandovici I et al. The imprinted Igf2-Igf2r axis is critical for matching placental microvasculature expansion to fetal growth. *Dev Cell.* (2022)
- Lopez-Tello, J et al. Ablation of PI3K-p110alpha Impairs Maternal Metabolic Adaptations to Pregnancy. *Frontiers in Cell and Developmental Biology* (2022)
- Lopez-Tello, J et al. Identification of Structural and Molecular Signatures Mediating Adaptive Changes in the Mouse Kidney in Response to Pregnancy. *International Journal of Molecular Sciences* (2022)
- Krilis M et al. Clinical relevance of nitrated beta 2-glycoprotein I in antiphospholipid syndrome: Implications for thrombosis risk. *J Autoimmun.* (2021)
- Mátyás, A et al. Hippocampal Sclerosis in Pilocarpine Epilepsy: Survival of Peptide-Containing Neurons and Learning and Memory Disturbances in the Adult NMRI Strain Mouse. *International Journal of Molecular Sciences* (2021)
- Trujillo-Cenoz O et al. The ependymal cell cytoskeleton in the normal and injured spinal cord of mice. *J Neurosci Res.* (2021)
- De Clercq K. et al. Double-label immunohistochemistry to assess labyrinth structure of the mouse placenta with stereology. *Placenta* (2020)
- Napso T, Hung YP, Davidge ST, Care AS, Sferruzzi-Perri AN. Advanced maternal age compromises fetal growth and induces sex-specific changes in placental phenotype in rats. *Sci Rep.* (2020)
- Aykroyd BRL, Tunster SJ, Sferruzzi-Perri AN. Igf2 deletion alters mouse placenta endocrine capacity in a sexually dimorphic manner. *J Endocrinol.* (2020)
- López-Tello J, Pérez-García V, Khaira J, et al. Fetal and trophoblast PI3K p110 α have distinct roles in regulating resource supply to the growing fetus in mice. *Elife.* (2019)
- Sozen B, Cox AL, De Jonghe J, et al. Self-Organization of Mouse Stem Cells into an Extended Potential Blastoid. *Dev Cell.* (2019)
- Kühnle A et al. Polysialic Acid Modulates the Binding of External Lactoferrin in Neutrophil Extracellular Traps. *Biology (Basel).* (2019)
- Huang MS, Lin WC, Chang JH, Cheng CH, Wang HY, Mou KY. The cysteine-free single mutant C32S of APEX2 is a highly expressed and active fusion tag for proximity labeling applications. *Protein Sci.* (2019)

- Napso, T et al. Advanced maternal age compromises fetal growth and induces sex-specific changes in placental phenotype in rats. *Scientific Reports* (2019)
- Elisa Balmas et al. Maternal group 2 innate lymphoid cells control fetal growth and protect from endotoxin-induced abortion in mice. *bioRxiv* (2018)
- Batareseh YS et al. Oleocanthol ameliorates amyloid- β oligomers' toxicity on astrocytes and neuronal cells: In vitro studies. *Neuroscience* (2017)
- Tatsukawa, H et al. Global identification and analysis of isozyme-specific possible substrates crosslinked by transglutaminases using substrate peptides in mouse liver fibrosis. *Scientific Reports* (2017)
- Vasek MJ, Garber C, Dorsey D, et al. A complement-microglial axis drives synapse loss during virus-induced memory impairment. *Nature*. (2016)
- Phillips et al. Entry Sites of Venezuelan and Western Equine Encephalitis Viruses in the Mouse Central Nervous System following Peripheral Infection. *Journal of Virology* (2016)
- Cantoni C, Bollman B, Licastro D, et al. TREM2 regulates microglial cell activation in response to demyelination in vivo. *Acta Neuropathol.* (2015)
- Matveev SV et al. A distinct subfraction of A β is responsible for the high-affinity Pittsburgh compound B-binding site in Alzheimer's disease brain. *Journal of Neurochemistry* (2014)
- LeVine H et al. Dihydroxybenzoic acid isomers differentially dissociate soluble biotinyl-A β (1-42) oligomers. *Biochemistry* (2012)
- Cowles CE, Li Y, Semmelhack MF, Cristea IM, Silhavy TJ. The free and bound forms of Lpp occupy distinct subcellular locations in *Escherichia coli*. *Mol Microbiol.* (2011)
- [View More ...](#)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.