

Datasheet for TW0020

Tween-20 Ultra Pure

Overview

Description:	Tween-20 Ultra Pure - TW0020
Item No.:	TW0020
Size:	50 mL

Product Details

Synonyms:	detergent
Detection Kit Type:	Colorimetric Western Blot Kit

Target Details

Relevant Links:	Tween 20 SDS
-----------------	--

Formulation

Physical State:	Liquid
Concentration:	neat

Shipping & Handling

Shipping Condition:	Ambient
Storage Condition:	Store container at room temperature (18° to 26° C) prior to opening.
Expiration:	Expiration date is two (2) years from date of receipt.

Images

**Bottle**

In biological techniques and sciences, Polysorbate 20 has a broad range of applications. One example, is its use as a washing agent in immunoassays, such as Western blots and ELISAs. It helps to prevent non-specific antibody binding. These buffers are used for washes between each immunoreactions, to remove unbound immunologicals, and eventually for incubation solutions of immunoreagents (labeled antibodies) to reduce unspecific background.

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

