

Datasheet for WIB-4625-005

Western Incubation Box

Overview

Description:	Western Blot Incubation Box, Large (4 5/8" X 3 1/2" X 1 1/8") , 5-PACK - WIB-4625-005
Item No.:	WIB-4625-005
Size:	5 Each
Applications:	WB

Product Details

Background:	Black incubation box designed for Western blotting, perfectly suited for use with fluorescence-conjugated antibodies, including DyLight conjugated secondary antibodies and IRDye antibody conjugates. This large box is best for volumes between 10 mL and 30 mL.
Synonyms:	Western Blot Box, Blot Box, incubation box

Application Details

Tested Applications:	WB
WB:	User Optimized

Formulation

Concentration:	5 Boxes
-----------------------	---------

Shipping & Handling

Shipping Condition:	Ambient
----------------------------	---------

Images

**Western Incubation Box**

Western Blot Incubation Box, Large (4 5/8" X 3 1/2" X 1 1/8"), 5-PACK

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.

