

LAB FAVORITES

VISIT WWW.ROCKLAND.COM FOR MORE

For over 60 years, Rockland has produced and manufactured antibodies and reagents in our USA-based facilities. Our products go through rigorous quality-control testing and multi-assay validation for consistent performance, with thousands of products in stock and ready to ship to your lab today.

Primary Antibodies

Product	Reactivity	Applications*	Item No.
AKT Antibody	Human, Mouse, Rat	IF, IHC, WB	100-401-401
ATM phospho S1981 Antibody	Human, Mouse, Rat	ELISA, FC, IF, WB	200-301-400
Beta Amyloid Antibody	Mouse	ELISA, IF, IHC, WB	600-401-253
Collagen Type I Antibody	Bovine, Human, Mouse, Rat, Pig	Dot Blot, IHC, WB	600-401-103-0.5
Collagen Type I Antibody Biotin Conjugated	Human, Mouse, Rat	Dot Blot, FC, Multiplex	600-406-103
Collagen Type II Antibody	Human	Dot Blot, IHC	600-401-104-0.5
Collagen Type III Antibody	Bovine, Human, Mouse, Rat, Pig	Dot Blot, IHC, WB	600-401-105-0.5
Collagen Type IV Antibody	Human	Dot Blot, IHC	600-401-106-0.5
Collagen Type VI Antibody	Human	Dot Blot, IHC	600-401-108-0.5
Human Red Blood Cell Antibody	Human	Agglutination, IHC	209-4139
IL-1 Beta Antibody	Human	WB	209-401-301
IL-6 Antibody	Human	Dot Blot, WB, IHC	109-401-310
Luciferase Antibody	Human	WB	200-101-C44
Myosin phospho S19/phospho S20 Antibody	Human	ELISA, IHC, WB	600-401-416
NOTCH 1 Antibody	Human	ELISA, IHC	100-401-407
Osteopontin Antibody	Human	ELISA, IHC, WB	100-401-404
Ovalbumin Antibody	Human	WB	200-401-033
Ovalbumin Antibody Peroxidase Conjugated	Rabbit	ELISA, WB	200-4333-0100
PPAR alpha Antibody	Mouse	ELISA, IF, IHC, WB	600-401-421
Rhodamine Antibody	Rhodamine	Dot Blot, ELISA, WB	200-301-246
ROBO-1 Antibody	Human	ELISA, IHC, WB	600-401-692
TNF alpha Antibody	Human, Mouse, Rat	IF, IHC, WB	209-401-306
Trypsin Inhibitor Antibody	Human	WB	200-401-079

*Applications shown include both tested and suggested. Please visit the product detail page on rockland.com for more information.

Epitope Tags

Product	Reactivity	Applications*	Item No.
Antibody for the detection of DYKDDDDK (FLAG™)	FLAG-Tag (C-Term)	WB, ELISA, IF, FC, ChIP, IP	200-301-383
Antibody for the detection of DYKDDDDK (FLAG™)	FLAG-Tag (C-Term)	WB, ELISA, IHC, IF, IP	600-401-383
Beta Galactosidase Antibody	b-GAL	WB, ELISA, IHC, IF	200-4136-0100
6X His Epitope Tag Antibody	His-Tag	ELISA, WB	200-301-382
6X His Tag Antibody HRP Conjugated	His-Tag	ELISA, WB	200-303-382
GFP Antibody	eGFP, rGFP	WB, ELISA, IHC, IF, FC, IP	600-101-215
GFP Antibody	eGFP, rGFP, YFP	WB, ELISA, IHC, IF, EM, IP	600-401-215
GFP Monoclonal Antibody	eGFP, rGFP	WB, ELISA, IHC, IF, FC, IP, ChIP	600-301-215
GST Antibody	GST-Tag	ELISA, WB	600-101-200

GST Antibody Biotin Conjugated	GST-Tag	ELISA	600-106-200
HA Epitope Tag Antibody	Human	ELISA, IHC, WB	600-401-384
RFP Antibody	RFP, rRFP	WB, IHC, IF, FC, EM	200-101-379
RFP Antibody	RFP, rRFP, tdTomato	WB, ELISA, IHC, IF, FC	600-901-379
RFP Antibody Pre-adsorbed	RFP, rRFP	WB, ELISA, IHC, IF, FC	600-401-379
RFP Antibody Biotin Conjugated Pre-adsorbed	RFP, mScarlet, rRFP, tdTomato	WB, ELISA, IHC, IF, FC	600-406-379

*Applications shown include both tested and suggested. Please visit the product detail page on rockland.com for more information.

Secondary Antibodies

Product	Host	Applications*	Item No.
Goat IgG (H&L) Antibody Peroxidase Conjugated	Donkey	WB, ELISA, IHC	605-703-002
Human IgG (H&L) Antibody Peroxidase Conjugated	Goat	WB, ELISA, Dot Blot	609-1302
Human IgG (H&L) Antibody Pre-Adsorbed	Goat	WB, ELISA, Dot Blot	609-101-123
Human IgG F(c) Antibody	Goat	WB, ELISA, LFA	109-1103
Mouse IgG (H&L) Antibody	Goat	WB, ELISA, LFA	610-1102
Mouse IgG (H&L) Antibody Peroxidase Conjugated	Goat	WB, ELISA, IHC, LFA	610-1302
Mouse IgG (H&L) Antibody Pre-adsorbed	Goat	WB, ELISA, LFA	610-101-121
Mouse IgG (H&L) Antibody	Rabbit	WB, ELISA, EM, IP	610-4102
Mouse IgG (H&L) Antibody Peroxidase Conjugated	Rabbit	ELISA, WB	610-4302
Rabbit IgG (H&L) Antibody	Goat	WB, ELISA, FC, Dot Blot, LFA	611-1102
Rabbit IgG (H&L) Antibody Alkaline Phosphatase Conjugated	Goat	WB, ELISA, IHC, Dot Blot	611-1502
Rabbit IgG (H&L) Antibody Pre-adsorbed	Goat	WB, ELISA, IF, Biochemical Assay, LFA	611-101-122
Rabbit IgG (H&L) Antibody CY5 Conjugated Pre-adsorbed	Goat	WB, ELISA, IF, Dot Blot, Microarray	611-110-122
Rabbit IgG (H&L) Antibody Peroxidase Conjugated	Goat	WB, ELISA, IHC	611-1302
Rabbit IgG (H&L) Antibody Peroxidase Conjugated Pre-adsorbed	Goat	WB, ELISA, IHC, Dot Blot	611-103-122
Rabbit IgG F(c) Antibody Pre-Adsorbed	Goat	WB, ELISA, Dot Blot	611-1123
Sheep IgG (H&L) Antibody	Rabbit	ELISA, LFA	613-4102

*Applications shown include both tested and suggested. Please visit the product detail page on rockland.com for more information.

TrueBlot® for IP/Western Blot

Over 400+ references! Decrease non-specific binding, background, and heavy/light chain interference for publication-quality images.

Product	Applications	Item No.
Mouse TrueBlot® ULTRA: Anti-Mouse Ig HRP	IP, WB	18-8817-33
Rabbit TrueBlot®: Anti-Rabbit IgG HRP	IP, WB	18-8816-31
TrueBlot® Anti-Rabbit Ig IP Agarose Bead	IP, WB	00-8800-25

Proteins & Peptides

Product	Applications	Item No.
Human Collagen Type I	WB	009-001-103
Human Collagen Type III	WB	009-001-105
Human Collagen Type VI	SDS-PAGE	009-001-108
Human IgG	ELISA, SDS-PAGE	009-0102
Human IgG F(ab') ₂	ELISA, SDS-PAGE	009-0104
Human Transferrin	SDS-PAGE	009-0134

Supporting Reagents

Product	Size(s)	Item No.
10X TTBS pH 7.5	50 mL, 1 L	MB-013
Blocking Buffer for Fluorescent Western Blotting	50 mL, 500 mL	MB-070
Bovine Serum Albumin - Fraction V	1 kg, 10 g, 50 g	BSA-50
Bovine Serum Albumin 30% Solution	10 mL, 50 mL, 500 mL, 1 L	BSA-30
Chemiluminescent FemtoMax™ Super Sensitive HRP Substrate	20 mL, 110 mL, 1 L	FEMTOMAX-110
ELISA Microwell Blocking Buffer with Stabilizer (Azide and Mercury Free)	100 mL, 1 L	MB-064-1000
Molecular Biology Grade UltraPure Water	10 x 100mL, 1 x 1 L	MB-010-1000
STREPTAVIDIN	SDS-PAGE	S000-01
STREPTAVIDIN ATTO 532 Conjugated	5 mg	S000-53
STREPTAVIDIN PEROXIDASE Conjugated	1 mg	S000-03
TMB ELISA PEROXIDASE SUBSTRATE	100 mL, 1 L	TMBE-1000
UltraPure Sterile Water	1 L, 4 x 1 L	MB-009-1000

Blood Products

Product	Size(s)	Item No.
Goat Serum (sterile)	50 mL, 100 mL, 250 mL, 500 mL, 1 L	D104-00-0050
Goat Serum Non-Sterile	50 mL, 100 mL, 250 mL, 500 mL, 1 L	D204-00-1000
Donkey Serum Non-Sterile	50 mL, 100 mL, 250 mL, 500 mL, 1 L	D220-00-0050



+1 484.791.3823
sales@rockland.com
www.rockland.com

